

Product datasheet for RC220903L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ApoER2 (LRP8) (NM_017522) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: ApoER2 (LRP8) (NM 017522) Human Tagged ORF Clone Lentiviral Particle

Symbol: ApoER2

Synonyms: APOER2; HSZ75190; LRP-8; MCI1

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_017522 **ORF Size:** 2100 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(RC220903).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 017522.4

 RefSeq Size:
 6994 bp

 RefSeq ORF:
 2103 bp

 Locus ID:
 7804

 UniProt ID:
 Q14114

 Cytogenetics:
 1p32.3

Domains: Idl_recept_b, EGF_CA, Idl_recept_a, EGF

Protein Families: Druggable Genome, Secreted Protein, Transmembrane






MW: 78.3 kDa

Gene Summary:

This gene encodes a member of the low density lipoprotein receptor (LDLR) family. Low density lipoprotein receptors are cell surface proteins that play roles in both signal transduction and receptor-mediated endocytosis of specific ligands for lysosomal degradation. The encoded protein plays a critical role in the migration of neurons during development by mediating Reelin signaling, and also functions as a receptor for the cholesterol transport protein apolipoprotein E. Expression of this gene may be a marker for major depressive disorder. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2011]

Product images:



[RC220903L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC220903L4V particle to overexpress human LRP8-mGFP fusion protein.