

## Product datasheet for **RC220908L4V**

### ATP12A (NM\_001676) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ATP12A (NM_001676) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ATP12A
Synonyms:	ATP1AL1; H-K-ATPase; HK
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001676
ORF Size:	3135 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220908).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001676.3</a> , <a href="#">NP_001667.3</a>
RefSeq Size:	3594 bp
RefSeq ORF:	3120 bp
Locus ID:	479
UniProt ID:	<a href="#">P54707</a>
Cytogenetics:	13q12.1-q12.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Oxidative phosphorylation



[View online »](#)

**MW:** 115.3 kDa

**Gene Summary:** The protein encoded by this gene belongs to the family of P-type cation transport ATPases. This gene encodes a catalytic subunit of the ouabain-sensitive H<sup>+</sup>/K<sup>+</sup> -ATPase that catalyzes the hydrolysis of ATP coupled with the exchange of H<sup>(+)</sup> and K<sup>(+)</sup> ions across the plasma membrane. It is also responsible for potassium absorption in various tissues. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]