

## Product datasheet for **RC221015**

### Dystrobrevin alpha (DTNA) (NM\_032975) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dystrobrevin alpha (DTNA) (NM_032975) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DTNA
Synonyms:	D18S892E; DRP3; DTN; DTN-A; LVNC1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC221015 representing NM\_032975  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGATTGAAGATAGTGGGAAAAGAGGAAATACCATGGCAGAAAGAAGACAGCTGTTTGCAGAGATGAGGG  
 CTAAGATCTGGATCGCATCCGACTCTCCACCTACAGAACAGCATGCAAGCTTAGGTTTGTTCAGAAGAA  
 ATGCAATTTGCACCTGGTGGACATATGGAATGTATAGAAGCATTGCGGGAAAATGCTCTGAACAACCTG  
 GACCCAAACACTGAACTCAACGTGTCCCCTTAGAGGCTGTGCTCTCCACTATTTTTTACCAGCTCAACA  
 AACGGATGCCAACACTCACAAATCCATGTGGAGCAGTCCATCAGCCTCCTCCTAACTCCTGCTTGC  
 AGCGTTTGTATCCGGAAGGCCATGGTAAATTTCAAGTATTTGCTGTCAAAATGGCTTTAGCCACATTTGT  
 GGAGGGAAGATCATGGACAAATTAAGATATATTTTCTCAATGATTTCTGACTCCAGTGGGTGATGGTTT  
 ATGGACGATATGACCAATTCCTTCGGGAAGTTCTCAAACACCCAGGCAGTTTTTGAAGGTCCTTCATT  
 TGTTTACACAGAACAGTCAGCCAGATCCTGTTTCTCCCAACAGAAAAAGTCACGTTAAATGGTTTCTTG  
 GACACGCTTATGTGAGATCCTCCCCCGAGTGTCTGGTCTGGTTGCCTCTCTGATCGACTAGCAAATG  
 TGGAAAATGTCTTCCATCCGGTTGAGTGTCTCTACTGCCACAGTGAAGTATGATGGGATTTTCGCTACCG  
 ATGCCAACAGTGTCAAAATACCAGCTCTGTCCAGGACTGCTTCTGGAGGGGACATGCCGGTGGTTTTCAT  
 AGCAACCAGCACCAATGAAAGAGTACACGTATGGAATCACCTGCTAAGAAGCTGACTAATGCATTA  
 GCAAGTCCCTGAGCTGTGCTTCCAGCCGTGAACCTTTCACCCCATGTTCCAGATCAGCCTGAGAAGCC  
 ACTCAACTGGCTCACATCGTTGATACTGGCCTCCCAGACCTGTAACCAGCATGAACGACACCCCTGTT  
 TCCCACTCTGTTCCCTCCTCAGGAAGTCCTTTTATTACCAGGAGCATGCTTGAGAGTTCAAACCGGCTTG  
 ATGAAGAACACAGGCTAATTGCCAGGTATGCGGCAAGGCTGGCAGCAGAGTCTCTTCGTCTCAGCCACC  
 TCAGCAGAGAAGTCTCCTGACATCTTTTACCACATCGATGCGAATAAGCAGCAAAGGCAGCTGATTGCT  
 GAGCTAGAAAACAAGAACAGAGAAATCTTACAGGAGATCCAGAGACTTCGGCTAGAGCATGAACAAGCTT  
 CTCAGCCACGCCAGAGAAGGCACAGCAAAACCCACCCCTGCTGGCAGAACTCCGGCTCCTCAGACAGCG  
 CAAAGATGAGCTGGAACAGAGAAATGTCTGCTCTCCAGGAGAGCCGGAGAGAGCTAATGGTCCAGTTGGAG  
 GGTCTCATGAAGCTACTAAAGACTCAGGGGGCAGGCTCTCCCGCTCCTCCCCAGCCACACCATCAGCA  
 GGCCAATTCCCATGCCATCCGGTCAGCGTCAGCCTGCTCCACCCGACGCACACGCCGAGGACTCCCT  
 CACAGGAGTAGGGGGAGATGTACAAGAGGCATTTGCACAAAGTTCAAGAAGAACTTAAGGAATGACTTG  
 CTAGTGGCTGCAGATTCATCACTAACACTATGTCTCTCTTGTGAAAGAGCTGAATTCGAGGTTGGGA  
 GTGAAACAGAGAGTAATGTGGATTCTGAATTTGCACGGACTCAGTTTGAAGATCTTGTCCCTCACCAAC  
 CTCTGAAAAGGCTTTTCTAGCGCAAATCCATGCCGAAAACCTGGGTACATTACAGTGGAGCTACCACA  
 AGTACCATGCGTGGCGACATGGTTACGGAGGATGCAGATCCCTATGTGCAGCCTGAAGATGAAAACATG  
 AAAATGACTCTGTCCGGCAGCTGGAGAATGAGCTCCAGATGGAGGAATACCTGAAACAGAAGCTGCAAGA  
 TGAAGCTTATCAGGTCAGCTTGAAGGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC221015 representing NM\_032975  
 Red=Cloning site Green=Tags(s)

```

MIEDSGKRGNTMAERRQLFAEMRAQDLDRIRLSTYRTACKLRFVQKKNLHLVDIWNVIEALRENALNNL
DPNTELNVSRLAVLSTIFYQLNKRMPPTHQIHVEQSISSLLNFLAADFPEGHGKISVFAVKMALATLC
GGKIMDKLRYIFSMISDSSGVMVYGRYDQFLREVLKLPVAVFEGPSFGYTEQSARSCFSQQKQVTLNGFL
DTLMSDPPPQCLVWLPLLHRLANVENVFHPVECSYCHSESMGFRYRCQQCHNYQLCQDCFWRGHAGGSH
SNQHQMKEYTSWKSPAKKLTNALSKSLSCASSREPLHPMFDPQPEKPLNLAHIVDTWPPRPVTSMNDTLF
SHSVPSGSPFITRSMLESSNRLDEEHRLIARYAARLAAESSSSQPPQQRSAVDISFTIDANKQQRLIA
ELENKNREILQEIQRLREHEQASQPTPEKAQQNPTLLAELRLLRQRKDELEQRMSALQESRRELMVQLE
GLMKLLKTQGAGSPRSSPHTISRPIMPIRSASACSTPTHTPQDSLTVGGDVQEAFQSSRRNLRNDL
LVAADSIINTMSSLVKELNSEVGETESNVDSEFARTQFEDLVPSPTSEKAFLAQIHARKPGYIHSGATT
STMGRGDMVTEDADPYVQPEDENYENDSVRQLENELQMEEYLKQKLQDEAYQVSLQG
  
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_032975

**ORF Size:** 2058 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_032975.3](#), [NP\\_116757.2](#)

**RefSeq Size:** 2511 bp

RefSeq ORF: 2061 bp

Locus ID: 1837

UniProt ID: [Q9Y4J8](#)

Cytogenetics: 18q12.1

Domains: ZnF\_ZZ

Protein Families: Druggable Genome

MW: 77.5 kDa

**Gene Summary:** The protein encoded by this gene belongs to the dystrobrevin subfamily of the dystrophin family. This protein is a component of the dystrophin-associated protein complex (DPC), which consists of dystrophin and several integral and peripheral membrane proteins, including dystroglycans, sarcoglycans, syntrophins and alpha- and beta-dystrobrevin. The DPC localizes to the sarcolemma and its disruption is associated with various forms of muscular dystrophy. Mutations in this gene are associated with left ventricular noncompaction with congenital heart defects. Multiple alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]