

## Product datasheet for **RC221032A1V**

### Human **DEFB108B (NM\_001002035) AAV Particle**

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human DEFB108B (NM\_001002035) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** DEFB108B  
**Synonyms:** DEFB-8; hBD-8  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC221032 representing NM\_001002035  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

**ATGAGGATTGCTGCTCCTCTTCGCCATTTCTTCTTTATGAGCCAAGTTCTACCAGCCAGGGGCAAAT**  
**TCAAGGAGATCTGTGAACGTCCAATGGCTCCTGTGCGGACTTTGCCTGAAACAGAAATCCATGTTGG**  
**GAGATGTTTAAATAGCCAACCTGCTGCCTGCCTCTGGGGCATCAACCAAGAATTGAGAGCACTACCCC**  
**AAAAAGGAC**

**ACGCGTACGCGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT**  
**ACAAGGATGACGACGATAAGGTTTAA**

**Protein Sequence:** >RC221032 representing NM\_001002035  
**Red**=Cloning site **Green**=Tags(s)

MRIAVLLFAIFFMSQVLPARGKFKEICERPNGSCRFLETEIHVGRCLNSQPCLPLGHQPRIESTTP  
KKD

**TRTRPLEQKLI SEEDLAANDILDYKDDDDKV**

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001002035  
**ORF Size:** 219 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001002035.1</u>
<b>RefSeq Size:</b>	222 bp
<b>RefSeq ORF:</b>	222 bp
<b>Locus ID:</b>	245911
<b>UniProt ID:</b>	<u>Q8NET1</u>
<b>Cytogenetics:</b>	11q13.4
<b>MW:</b>	8.1 kDa