

## Product datasheet for **RC221045L1V**

### ADAM21 (NM\_003813) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAM21 (NM_003813) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAM21
Synonyms:	ADAM 21; ADAM31
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003813
ORF Size:	2166 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221045).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003813.2</a>
RefSeq Size:	2406 bp
RefSeq ORF:	2169 bp
Locus ID:	8747
UniProt ID:	<a href="#">Q9UKJ8</a>
Cytogenetics:	14q24.2
Protein Families:	Protease, Transmembrane
MW:	80.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The expression of this gene expression is testis-specific. [provided by RefSeq, May 2011]