

Product datasheet for **RC221048L1V**

PPP2R2D (NM_001003656) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP2R2D (NM_001003656) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPP2R2D
Synonyms:	MDS026
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001003656
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221048).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001003656.1 , NP_001003656.1
RefSeq Size:	2091 bp
RefSeq ORF:	680 bp
Locus ID:	55844
Cytogenetics:	10q26.3
Protein Families:	Druggable Genome, Phosphatase
Protein Pathways:	Tight junction
MW:	26.1 kDa



[View online »](#)

Gene Summary:

B regulatory subunit of protein phosphatase 2A (PP2A) that plays a key role in cell cycle by controlling mitosis entry and exit. The activity of PP2A complexes containing PPP2R2D (PR55-delta) fluctuate during the cell cycle: the activity is high in interphase and low in mitosis. During mitosis, activity of PP2A is inhibited via interaction with phosphorylated ENSA and ARPP19 inhibitors. Within the PP2A complexes, the B regulatory subunits modulate substrate selectivity and catalytic activity, and also may direct the localization of the catalytic enzyme to a particular subcellular compartment (By similarity).[UniProtKB/Swiss-Prot Function]