

## Product datasheet for **RC221052L3V**

### Kallikrein 8 (KLK8) (NM\_144506) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Kallikrein 8 (KLK8) (NM_144506) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Kallikrein 8
Synonyms:	HNP; NP; NRPN; PRSS19; TADG14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144506
ORF Size:	357 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221052).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_144506.1</a>
RefSeq Size:	590 bp
RefSeq ORF:	360 bp
Locus ID:	11202
UniProt ID:	<a href="#">O60259</a>
Cytogenetics:	19q13.41
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	12.5 kDa



[View online »](#)

**Gene Summary:**

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in tandem in a gene cluster on chromosome 19. The encoded protein may be involved in proteolytic cascade in the skin and may serve as a biomarker for ovarian cancer. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2013]