

Product datasheet for **RC221076**

OR5D13 (NM_001001967) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR5D13 (NM_001001967) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR5D13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC221076 representing NM_001001967 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGCATCTGAAAGAAATCAAAGCAGCACACCCACTTTTATTCTCTTGGGTTTTTCAGAATACCCAG
AAATCCAGGTTCCACTCTTTCTGGTTTTCTTGTCGTCTACACAGTCACTGTAGTGGGAACTTGGGCAT
GATAATAATCATCAGACTCAATTCAAAACCTCCATACAATCATGTGCTTTTCTTAGTCACTTGTCTTG
ACAGACTTCTGTTTTCCACTGTAGTTACACCTAACTGTTGGAGAACTTGGTTGTGGAATACAGAACCA
TCTTTTCTCTGGTTGCATCATGCAATTTGTTTTGCTTGCAATTTTGGAGTGACAGAACTTTTCATGTT
AGCAGCGATGGCTTATGACCGTTTTGTGGCAGTTTGTAACCCCTTGCTGTATACCACTATTATGTCTCAG
AAGCTCTGTGCTCTTCTGGTGGCTGGGTCTATACATGGGGGATAGTGTGCTCCCTGATACTCACATATT
TTCTTCTTGACTTATCGTTTTGTGAATCTACCTTCATAAATAATTTTATCTGTGACCACTCTGTAATTGT
TTCTGCCTCCTACTCAGACCCCTATATCAGCCAGAGGCTATGCTTTATTATTGCCATATTCAATGAGGTG
AGCAGCCTAATTATCATTCTGACATCATATATGCTTATTTTCACTACCATTATGAAGATGCGATCTGCAA
GTGGGCGCCAGAAAACCTTCTCCACCTGTGCCTCCACCTGACAGCCATCACTATCTTCCATGGAACAT
CCTTTTCCTTTACTGTGTTCTAATCCTAAAACCTCTAGCCTCATAGTTACAGTGGCTTCTGTGTTTTAC
ACAGTGGCGATTCCAATGCTGAACCCATTGATCTACAGCCTTAGGAACAAAGATATCAATAACATGTTTG
AAAAATTAGTTGTCACCAAATTGATTTACCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC221076 representing NM_001001967
 Red=Cloning site Green=Tags(s)

MMASERNQSSTPTFILLGFSEYPEIQVPLFLVFLFVYTVTVVGNLGMIIIRLNSKLHTIMCFFLSHLSL
 TDFCFSTVVTPKLEENLVVEYRTISFSGCIMQFCFACIFGVTETFMFLAAMAYDRFVAVCKPLLYTTIMSQ
 KLCALLVAGSYTWGIVCSLILTYFLDLFSCESTFINNFICDHSVIVSASYSDPYISQRLCFIIAIFNEV
 SSLIIILT SYMLIFTTIMKMRASGRQKTFSTCASHLTAITIFHGTLILFLYCVPNPKTSSLIVTVASVFY
 TVAIPMLNPLIYSLRNKDINNMFEEKLVVTKLIYH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8001_e09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001001967

ORF Size: 942 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001001967.1, NP_001001967.1</u>
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	390142
UniProt ID:	<u>Q8NGL4</u>
Cytogenetics:	11q11
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	35.3 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. This olfactory receptor gene is a segregating pseudogene, where some individuals have an allele that encodes a functional olfactory receptor, while other individuals have an allele encoding a protein that is predicted to be non-functional. [provided by RefSeq, Jun 2015]

Product images:



Circular map for RC221076