

## Product datasheet for **RC221083L4V**

### **TBC1D10A (NM\_031937) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TBC1D10A (NM_031937) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TBC1D10A
Synonyms:	dj130H16.1; dj130H16.2; EPI64; TBC1D10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031937
ORF Size:	1524 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221083).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031937.1</a>
RefSeq Size:	1973 bp
RefSeq ORF:	1527 bp
Locus ID:	83874
UniProt ID:	<a href="#">Q9BXI6</a>
Cytogenetics:	22q12.2
MW:	56.9 kDa


[View online »](#)

**Gene Summary:**

Acts as GTPase-activating protein for RAB27A, but not for RAB2A, RAB3A, nor RAB4A.  
[UniProtKB/Swiss-Prot Function]