

## Product datasheet for **RC221086L3V**

### CCDC132 (VPS50) (NM\_024553) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC132 (VPS50) (NM_024553) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VPS50
Synonyms:	CCDC132; VPS54L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024553
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221086).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024553.2</a> , <a href="#">NP_078829.1</a>
RefSeq Size:	1216 bp
RefSeq ORF:	984 bp
Locus ID:	55610
UniProt ID:	<a href="#">Q96JG6</a>
Cytogenetics:	7q21.2-q21.3
MW:	37.3 kDa


[View online »](#)

**Gene Summary:**

Acts as component of the EARP complex that is involved in endocytic recycling. The EARP complex associates with Rab4-positive endosomes and promotes recycling of internalized transferrin receptor (TFRC) to the plasma membrane. Within the EARP complex, required to tether the complex to recycling endosomes. Not involved in retrograde transport from early and late endosomes to the trans-Golgi network (TGN).[UniProtKB/Swiss-Prot Function]