

Product datasheet for **RC221090A1V**

Human SNHG32 (NM_001040437) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNHG32 (NM_001040437) AAV Particle
Tag: Myc-DDK
Symbol: SNHG32
Synonyms: D6S57; G8
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC221090 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGAGAAGCTTTGTATGGCTGTCATGCTTAGACAGTGATTCTGCAACTTGACCTTCAGGCTGGGAG
AGGTGGAGAGCCATGCCTGTTCTCCTTCCTTGCATGGAATTTGCTGACACAATATCTCCGCCTGGTGC
TGGGCATATCCTAAGAATTACAACCTTCCTGTATTATCCTGTGTGAGCAGCTGCACCTATTGGGGGA
AAAATGCCTGAAAAT

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC221090 protein sequence
 Red=Cloning site Green=Tags(s)

MERSFVWL SCLDSDSNLTFRLGEVESHACSPSLLWNLLTQYLPPGAGHILRTYNFPVLSVSSCHLIGG
 KMPEN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001040437
ORF Size: 225 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001040437.1](#), [NP_001035527.1](#)

RefSeq Size: 845 bp

RefSeq ORF: 228 bp

Locus ID: 50854

Cytogenetics: 6p21.33

MW: 8.3 kDa