

## Product datasheet for **RC221125L1V**

### **MEX3C (NM\_016626) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MEX3C (NM_016626) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MEX3C
Synonyms:	BM-013; MEX-3C; RKHD2; RNF194
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016626
ORF Size:	1977 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221125).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016626.3</a>
RefSeq Size:	3797 bp
RefSeq ORF:	1980 bp
Locus ID:	51320
UniProt ID:	<a href="#">Q5U5Q3</a>
Cytogenetics:	18q21.2
Domains:	RING
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 69.2 kDa

**Gene Summary:** This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to hypertension. [provided by RefSeq, Oct 2009]