

## Product datasheet for **RC221139L4V**

### ADAM2 (NM\_001464) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAM2 (NM_001464) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAM2
Synonyms:	CRYN1; CRYN2; CT15; FTNB; PH-30b; PH30; PH30-beta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001464
ORF Size:	2205 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221139).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001464.3</a>
RefSeq Size:	2657 bp
RefSeq ORF:	2208 bp
Locus ID:	2515
UniProt ID:	<a href="#">Q99965</a>
Cytogenetics:	8p11.22
Protein Families:	Druggable Genome, Protease, Transmembrane
MW:	82.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The encoded protein is a subunit of an integral sperm membrane glycoprotein called fertilin, which plays an important role in sperm-egg interactions. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2013]