

Product datasheet for **RC221177L4V**

NICE2 (S100A7A) (NM_176823) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NICE2 (S100A7A) (NM_176823) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NICE2
Synonyms:	NICE-2; NICE2; S100A7f; S100A7L1; S100A15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_176823
ORF Size:	303 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_176823.2
RefSeq Size:	4351 bp
RefSeq ORF:	306 bp
Locus ID:	338324
UniProt ID:	Q86SG5
Cytogenetics:	1q21.3
MW:	11.1 kDa


[View online »](#)

Gene Summary:

May be involved in epidermal differentiation and inflammation and might therefore be important for the pathogenesis of psoriasis and other diseases.[UniProtKB/Swiss-Prot Function]