

## Product datasheet for **RC221210L4V**

### **GPR126 (ADGRG6) (NM\_198569) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GPR126 (ADGRG6) (NM_198569) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR126
Synonyms:	APG1; DREG; GPR126; LCCS9; PR126; PS1TP2; STQTL1; VIGR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198569
ORF Size:	3750 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221210).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198569.2</a> , <a href="#">NP_940971.1</a>
RefSeq Size:	6934 bp
RefSeq ORF:	3753 bp
Locus ID:	57211
UniProt ID:	<a href="#">Q86SQ4</a>
Cytogenetics:	6q24.2
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	139.9 kDa



[View online »](#)

**Gene Summary:**

This gene, which is upregulated in human umbilical vein endothelial cells, encodes a G protein-coupled receptor. Variations in this gene can affect a person's stature. Multiple transcript variants encoding different proteins have been found for this gene. [provided by RefSeq, Mar 2009]