

Product datasheet for **RC221278A1V**

Human LILRA5 (NM_181986) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LILRA5 (NM_181986) AAV Particle
Tag:	Myc-DDK
Symbol:	LILRA5
Synonyms:	CD85; CD85F; ILT-11; ILT11; LILRB7; LIR-9; LIR9
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221278 representing NM_181986 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCACCATGGTCTCATCCATCTGCACAGCTGCAGCCAGTGGGAGGAGACGCCGTGAGCCCTGCCCTCA
TGGTTCGTCTGCCTCGGGAACCTCTCAAAGCCACCCTCTGGGCTGAGCCAGGCTCTGTGATCAGCCG
GGGGAACCTCTGTGACCATCCGGTGTGAGGGGACCCTGGAGGCCAGGAATACCGTCTGGTTAAAGAGGGA
AGCCCAGAACCCTGGGACACACAGAACCCTGGAGCCCAAGAACAAGGCCAGATTCTCCATCCCATCCA
TGACAGAGCACCATGCAGGGAGATACCGCTGTTACTACTACAGCCCTGCAGGCTGGTCAGAGCCCAGCGA
CCCCCTGGAGCTGGTGGTGACAGGATTCTACAACAAACCCACCCTCTCAGCCCTGCCAGTCCTGTGGTG
ACCTCAGGAGAGAACGTGACCCTCCAGTGTGGCTCACGGCTGAGATTCGACAGGTTATTCTGACTGAGG
AAGGAGACCACAAGCTCTCCTGGACCTTGGAATCACAGCTGACCCCAAGTGGGAGTTCCAGGCCCTGTT
CCCTGTGGGCCCTGTGACCCCAAGCCAGGTGGATGCTCAGATGCTATGGCTCTCGCAGGCATATCCTG
CAGGTATGGTCAGAACCCAGTGACCTCCTGGAGATTCCGGTCTCAGGTGAGGAAGCCACAGTCTTCTCTA
GTACAATTCAGGAAGCCAGACAGGTTGTGGAGAGCTTTACAGGCAGGGCAGCCCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC221278 representing NM_181986 Red=Cloning site Green=Tags(s) MAPWSPSAQLQPVGDAVSPALMVLLCLGNLSKATLWAEPGSVISRGNSVTIRCQGTLEAQEYRLVKEG SPEPWTQNPLEPKNKARFSIPSMTEHHAGRYRCYYSPAGWSEPSDPLELVVTGFYNKPTLSALPSPVV TSGENVTLQCGSRLRFDRFILTEGDHKLWTLDSQLTPSGQFQALFPVGPVTPSHRWMLRCYGSRRHIL QVWSEPSDLLLEIPVSGEEATVFSSTIQGSQTGCGELYRQGSPC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_181986
ORF Size:	759 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_181986.1
RefSeq Size:	1004 bp
RefSeq ORF:	762 bp
Locus ID:	353514
UniProt ID:	A6NI73
Cytogenetics:	19q13.42
MW:	27.7 kDa