

Product datasheet for **RC221298L4V**

PACAP (ADCYAP1) (NM_001099733) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PACAP (ADCYAP1) (NM_001099733) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PACAP
Synonyms:	PACAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001099733
ORF Size:	528 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221298).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001099733.1
RefSeq Size:	3187 bp
RefSeq ORF:	531 bp
Locus ID:	116
UniProt ID:	P18509
Cytogenetics:	18p11.32
Protein Families:	Secreted Protein
MW:	18.8 kDa



[View online »](#)

Gene Summary:

This gene encodes a secreted proprotein that is further processed into multiple mature peptides. These peptides stimulate adenylate cyclase and increase cyclic adenosine monophosphate (cAMP) levels, resulting in the transcriptional activation of target genes. The products of this gene are key mediators of neuroendocrine stress responses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2013]