

## Product datasheet for **RC221301L3V**

### **NDRG4 (NM\_022910) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NDRG4 (NM_022910) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDRG4
Synonyms:	BDM1; SMAP-8; SMAP8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022910
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221301).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022910.1</a> , <a href="#">NP_075061.1</a>
RefSeq Size:	3453 bp
RefSeq ORF:	1116 bp
Locus ID:	65009
UniProt ID:	<a href="#">Q9ULP0</a>
Cytogenetics:	16q21
Domains:	Ndr
MW:	40.6 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that is required for cell cycle progression and survival in primary astrocytes and may be involved in the regulation of mitogenic signalling in vascular smooth muscles cells. Alternative splicing results in multiple transcripts encoding different isoforms.[provided by RefSeq, Jun 2011]