

Product datasheet for **RC221305L1V**

FER1L6 (NM_001039112) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FER1L6 (NM_001039112) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FER1L6
Synonyms:	C8ORFK23
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001039112
ORF Size:	5571 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221305).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039112.1 , NP_001034201.1
RefSeq Size:	6058 bp
RefSeq ORF:	5574 bp
Locus ID:	654463
UniProt ID:	Q2WGI9
Cytogenetics:	8q24.13
Protein Families:	Transmembrane
MW:	209.1 kDa



[View online »](#)