

Product datasheet for **RC221307L3V**

C22orf24 (NM_015372) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C22orf24 (NM_015372) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C22orf24
Synonyms:	HSN44A4A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015372
ORF Size:	480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015372.1 , NP_056187.1
RefSeq Size:	1247 bp
RefSeq ORF:	483 bp
Locus ID:	25775
UniProt ID:	Q9Y442
Cytogenetics:	22q12.3
Protein Families:	Transmembrane
MW:	17.6 kDa



[View online »](#)