

Product datasheet for **RC221309L1V**

PKC nu (PRKD3) (NM_005813) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PKC nu (PRKD3) (NM_005813) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PKC nu
Synonyms:	EPK2; nPKC-NU; PKC-NU; PKD3; PRKCN
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005813
ORF Size:	2670 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221309).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005813.3
RefSeq Size:	5907 bp
RefSeq ORF:	2673 bp
Locus ID:	23683
UniProt ID:	O94806
Cytogenetics:	2p22.2
Domains:	pkinaase, TyrKc, PH, DAG_PE-bind, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 100.3 kDa

Gene Summary: This gene belongs to the multigene protein kinase D family of serine/threonine kinases, which bind diacylglycerol and phorbol esters. Members of this family are characterized by an N-terminal regulatory domain comprised of a tandem repeat of cysteine-rich zinc-finger motifs and a pleckstrin domain. The C-terminal region contains the catalytic domain and is distantly related to calcium-regulated kinases. Catalytic activity of this enzyme promotes its nuclear localization. This protein has been implicated in a variety of functions including negative regulation of human airway epithelial barrier formation, growth regulation of breast and prostate cancer cells, and vesicle trafficking. [provided by RefSeq, Jan 2015]