

Product datasheet for **RC221326A1V**

Human Marapsin (PRSS27) (NM_031948) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Marapsin (PRSS27) (NM_031948) AAV Particle
Tag:	Myc-DDK
Symbol:	Marapsin
Synonyms:	CAPH2; MPN
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221326 representing NM_031948 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGGCGGCCGGCGCGGTGCCGCTCCTGCTGCTGTGTGTTTTGGGTCTCAGAGGGCCAAGGCAGCAA
CAGCCTGTGGTCGCCCCAGGATGCTGAACCGAATGGTGGGCGGGCAGGACACGCAGGAGGGCGAGTGGCC
CTGGCAAGTCAGCATCCAGCGCAACGGAAGCCACTTCTGCGGGGGCAGCCTCATCGCGGAGCAGTGGGTC
CTGACGGCTGCGCACTGCTTCCGCAACACCTCTGAGACGTCCTGTACCAGGTCCTGCTGGGGCAAGGC
AGCTAGTGCAGCCGGGACCACACGCTATGTATGCCCGGGTGAGGCAGGTGGAGAGCAACCCCTGTACCA
GGGCACGGCCTCCAGCGCTGACGTGGCCCTGGTGGAGCTGGAGGCACCAAGTGGCCCTTACCAATTACATC
CTCCCCGTGTGCCTGCCTGACCCCTCGGTGATCTTTGAGACGGGCATGAACTGCTGGGTCACTGGCTGGG
GCAGCCCCAGTGAGGAAGACCTCCTGCCGAACCGCGGATCCTGCAGAACTCGCTGTGCCCATCATCGA
CACACCAAGTGCAACCTGCTCTACAGCAAAGACACCGAGTTTGGCTACCAACCCAAAACCATCAAGAAT
GACATGCTGTGCGCCGGCTTCGAGGAGGGCAAGAAGGATGCCTGCAAGGGCGACTCGGGCGGCCCTGG
TGTGCCTCGTGGGTGAGTCGTGGCTGCAGGCGGGGTGATCAGCTGGGTGAGGGCTGTGCCCGCCAGAA
CCGCCCAGGTGTCTACATCCGTGTACCGCCACCACAACCTGGATCCATCGGATCATCCCCAACTGCAG
TTCCAGCCAGCGAGGTTGGGCGGCCAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC221326 representing NM_031948 Red=Cloning site Green=Tags(s) MRRPAAVPLLLLLCFGSQRAKAATACGRPRMLNRMVGGQDTQEGEWPWQVSIQRNGSHFCGGSLIAEQWV LTAAHCFRNTSETSLYQVLLGARQLVQPGPHAMYARVRQVESNPLYQGTASSADVALVELEAPVPFTNYI LPVCLPDPSVIFETGMNCWVTGWGSPSEEDLLPEPRILQKLAVPIIDTPKCNLLYSKDTEFGYQPKTIKN DMLCAGFEEGKKDACKGDSGGPLVCLVGQSWLQAGVISWGEGCARQNRPGVYIRVTAHHNWIHRIIPKLQ FQPARLGGQK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031948
ORF Size:	870 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031948.3</u>
RefSeq Size:	1462 bp
RefSeq ORF:	873 bp
Locus ID:	83886
UniProt ID:	<u>Q9BQR3</u>
Cytogenetics:	16p13.3
MW:	31.8 kDa