

Product datasheet for **RC221352L4V**

TAS2R42 (NM_181429) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TAS2R42 (NM_181429) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TAS2R42
Synonyms:	T2R24; T2R42; T2R55; TAS2R55
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181429
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221352).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181429.1 , NP_852094.1
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	353164
Cytogenetics:	12p13.2
Protein Families:	Transmembrane
Protein Pathways:	Taste transduction
MW:	36.2 kDa



[View online »](#)