

Product datasheet for **RC221395L4V**

TMPRSS11A (NM_182606) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMPRSS11A (NM_182606) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMPRSS11A
Synonyms:	ECRG1; HATL1; HESP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182606
ORF Size:	1263 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221395).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182606.2 , NP_872412.2
RefSeq Size:	2429 bp
RefSeq ORF:	1266 bp
Locus ID:	339967
UniProt ID:	Q6ZMR5
Cytogenetics:	4q13.2
Protein Families:	Druggable Genome, Transmembrane
MW:	47.4 kDa



[View online »](#)

Gene Summary:

Probable serine protease which may play a role in cellular senescence. Overexpression inhibits cell growth and induce G1 cell cycle arrest.[UniProtKB/Swiss-Prot Function]