

Product datasheet for **RC221411L4V**

SENP7 (NM_020654) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SENP7 (NM_020654) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SENP7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020654
ORF Size:	3150 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221411).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020654.3
RefSeq Size:	5018 bp
RefSeq ORF:	3153 bp
Locus ID:	57337
UniProt ID:	Q9BQF6
Cytogenetics:	3q12.3
Domains:	Peptidase_C48
Protein Families:	Druggable Genome, Protease
MW:	119.7 kDa



[View online »](#)

Gene Summary:

The reversible posttranslational modification of proteins by the addition of small ubiquitin-like SUMO proteins (see SUMO1; MIM 601912) is required for many cellular processes. SUMO-specific proteases, such as SENP7, process SUMO precursors to generate a C-terminal diglycine motif required for the conjugation reaction. They also display isopeptidase activity for deconjugation of SUMO-conjugated substrates (Lima and Reverter, 2008 [PubMed 18799455]).[supplied by OMIM, Jun 2009]