

Product datasheet for **RC221424L3V**

AMD1 (NM_001634) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AMD1 (NM_001634) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AMD1
Synonyms:	ADOMETDC; AMD; SAMDC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001634
ORF Size:	1002 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221424).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001634.4
RefSeq Size:	3421 bp
RefSeq ORF:	1005 bp
Locus ID:	262
UniProt ID:	P17707
Cytogenetics:	6q21
Domains:	SAM_decarbox
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Arginine and proline metabolism, Cysteine and methionine metabolism, Metabolic pathways

MW: 38.2 kDa

Gene Summary: This gene encodes an important intermediate enzyme in polyamine biosynthesis. The polyamines spermine, spermidine, and putrescine are low-molecular-weight aliphatic amines essential for cellular proliferation and tumor promotion. Multiple alternatively spliced transcript variants have been identified. Pseudogenes of this gene are found on chromosomes 5, 6, 10, X and Y. [provided by RefSeq, Dec 2013]