

Product datasheet for **RC221429A1V**

Human DEPDC4 (NM_152317) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DEPDC4 (NM_152317) AAV Particle
Tag:	Myc-DDK
Symbol:	DEPDC4
Synonyms:	DEP.4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221429 representing NM_152317 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGTGCCAGGGGAGGAGCCAGCGCGCAGCTTATGGCGGTTCTTTGACTCCGAGGTTCCGTAGACTTG
TCAGTCAGAACGAGCTTCCGGGCCAGGGCTGAACGGGCCAAGTTCAGGAACCGTAGAGATGGCTTCTG
CCGAAAAGGAGGACAGGATGCTCTGGTCCTTTCAAGCTACTCAGCTATGGGATGGTATTATCACTCT
CTTCAGGCCCAAGTGGAATAAAAAGAAGGAGGCATCATTTACAGACATACAAAGACTGTTTCACTGGTT
CTGATGCTGTCGATGTGGTCTTAAGTCATCTTATGCAAAACACGTGCCTAAGTAGCAATGACATCTTTG
TCTTAAAGGGGTTTCATCTTTGCCAAGTTCTAATGAATCACAAGTATTTGAACCAAGTAGGAATGAAGAAG
CTTTTCAAAAAGAAAAGGAATTAGAATTTGAAGATTCCAACATTAGTCTCTACAGGTTTCTAGGCAATA
AATCATCTTACGATTGTTGCAAAAGACAAAAGGATGCTGAGAATGAGTTCAATGAAACCTTGAGACCAGG
ATATGAAATGATTTCAAATCCTCTAGCACAGGAGATTGGTGAGGAAAGAATTGAGGAACCTTATTCATACA
ATAAATGGGAATCCAGCTTTATGTCCAATATCACAGTTCAGAAACCTTTCTCCGGCTTTCAAAAGAAG
ATGTTTGAAAGAACAACATTATTATGTCTTCTCAATTGATTCACCTCCATTCTGGACAATATTTT
GGAGCCTCCAGTTAAACACAAAATCTCAACTAAACAAAGAGGAAGATCTTGTATCACTAACACTTGC
CTAGACAGAGAACTTATTCCAAGCTTATGTCTACCTGAAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC221429 representing NM_152317 Red=Cloning site Green=Tags(s) MYPGEEPARELMAVLLTPRFRLVSQNELPGPGLNGPSSRNRRDGFCKRRTGCSGPFQATQLWDGIHS LQAQVEIKRRRHHLQTYKDCFTGSDAVDVVLSHLMQNTCLSSNDISCLKGVHLCQVLMNHKVFEPVGMKK LFKKEKELEFEDSNISLYRFLGNKSSYDCCRQKDAENEFNETLRPGYEMISNPLAQEIGEERIEELIHT INGNPALCPNITVQKPFLRLSKEDVWKEQTLLCLLQLIHLPLDNILEPPVKTQNLQLNKEEDLVITNTC LDRELIPSLCLPEK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152317
ORF Size:	882 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152317.1</u>
RefSeq Size:	2337 bp
RefSeq ORF:	885 bp
Locus ID:	120863
UniProt ID:	<u>Q8N2C3</u>
Cytogenetics:	12q23.1
MW:	33.5 kDa