

Product datasheet for **RC221436A1V**

Human Kallikrein 6 (KLK6) (NM_001012965) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Kallikrein 6 (KLK6) (NM_001012965) AAV Particle
Tag:	Myc-DDK
Symbol:	Kallikrein 6
Synonyms:	Bssp; hK6; Klk7; PRSS9; PRSS18; SP59
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221436 representing NM_001012965 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGAAGCTGATGGTGGTGTGAGTCTGATTGCTGCAGCCTGGGCAGAGGAGCAGAATAAGTTGGTGC
ATGGCGGACCCTGCGACAAGACATCTACCCCTACCAAGCTGCCCTCTACACCTCGGGCCACTTGCTCTG
TGGTGGGGTCTTATCCATCCACTGTGGGTCTCACAGCTGCCCACTGCAAAAAACCGAATCTTCAGGTC
TTCCTGGGGAAGCATAACCTTCGGCAAAGGGAGAGTTCACAGGAGCAGAGTCTGTTGTCCGGGCTGTGA
TCCACCCTGACTATGATGCCGCCAGCCATGACCAGGACATCATGCTGTTGCGCCTGGCAGCCAGCCAA
ACTCTCTGAACATCCAGCCCTTCCCCTGGAGAGGGACTGCTCAGCCAACACCACGCTGCCACATC
CTGGGCTGGGGCAAGACAGCAGATGGTATTTCCCTGACACCATCCAGTGTGCATACATCCACCTGGTGT
CCCGTGAGGAGTGTGAGCATGCCTACCCTGGCCAGATCACCCAGAACATGTTGTGTGCTGGGGATGAGAA
GTACGGGAAGGATTCTGCCAGGGTATTCTGGGGTCCGCTGGTATGTGGAGACCACCTCCGAGGCCTT
GTGTCATGGGGTAACATCCCCTGTGGATCAAAGGAGAAGCCAGGAGTCTACACCAACGTCTGCAGATACA
CGAACTGGATCCAAAAACCATTCAGGCCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC221436 representing NM_001012965 Red=Cloning site Green=Tags(s) MKKLMVVL SLIAAAWAEQNKLVHGGPCDKTSHPYQAALYTSGHLLCGGVLIHPLWVLTAAHCKKPNLQV FLGKHNL RQRESSQEQQSSVVRVHPDYDAASHDQDIMLLRLARPAKLSLIQPLPLERDCSANTTSCHI LGWGKTADGDFPDTIQCAYIHLVSREECEHAYPGQITQNMLCAGDEKYGKDSCQGDSSGGLVCGDHLRGL VSWGNI PCGSKEKPGVYTNVCRYTNWIQKTIQAK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001012965
ORF Size:	732 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001012965.1</u>
RefSeq Size:	1495 bp
RefSeq ORF:	414 bp
Locus ID:	5653
UniProt ID:	<u>Q92876</u>
Cytogenetics:	19q13.41
MW:	26.86 kDa