

Product datasheet for **RC221437L1V**


LRRC15 (NM_130830) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LRRC15 (NM_130830) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LRRC15
Synonyms:	LIB
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_130830
ORF Size:	1743 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221437).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130830.2
RefSeq Size:	5883 bp
RefSeq ORF:	1746 bp
Locus ID:	131578
UniProt ID:	Q8TF66
Cytogenetics:	3q29
Protein Families:	Transmembrane
MW:	64.4 kDa



[View online »](#)

Product images:

[RC221437L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC221437L1V particle to overexpress human LRRC15-Myc-DDK fusion protein.