

Product datasheet for RC221441L4

C15ORF27 (TMEM266) (NM_152335) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C15ORF27 (TMEM266) (NM_152335) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	C15ORF27
Synonyms:	C15orf27; HshVVRP1; hTMEM266; HVRP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221441).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_152335
ORF Size:	1593 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152335.2
RefSeq Size:	2414 bp
RefSeq ORF:	1596 bp
Locus ID:	123591
UniProt ID:	Q2M3C6
Cytogenetics:	15q24.2
Protein Families:	Druggable Genome, Transmembrane
MW:	58.3 kDa
Gene Summary:	Voltage-sensor protein present on the post-synaptic side of glutamatergic mossy fibers and granule cells in the cerebellum (PubMed:25165868, PubMed:30810529). Despite the presence of a voltage-sensor segment, does not form a functional ion channel and its precise role remains unclear (PubMed:25165868, PubMed:30810529). Undergoes both rapid and slow structural rearrangements in response to changes in voltage (PubMed:30810529). Contains a zinc-binding site that can regulate the slow conformational transition (PubMed:30810529). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC221441L4