

Product datasheet for **RC221450L3V**

NSFL1C (NM_018839) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NSFL1C (NM_018839) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NSFL1C
Synonyms:	dj776F14.1; P47; UBX1; UBXD10; UBXN2C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018839
ORF Size:	1017 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221450).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018839.3 , NP_061327.2
RefSeq Size:	3551 bp
RefSeq ORF:	1020 bp
Locus ID:	55968
UniProt ID:	Q9UNZ2
Cytogenetics:	20p13
Domains:	UBX, FAF
MW:	37.1 kDa



[View online »](#)

Gene Summary:

N-ethylmaleimide-sensitive factor (NSF) and valosin-containing protein (p97) are two ATPases known to be involved in transport vesicle/target membrane fusion and fusions between membrane compartments. A trimer of the protein encoded by this gene binds a hexamer of cytosolic p97 and is required for p97-mediated regrowth of Golgi cisternae from mitotic Golgi fragments. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8. [provided by RefSeq, May 2011]