

Product datasheet for **RC221458L3V**

UHRF1BP1 (NM_017754) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UHRF1BP1 (NM_017754) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UHRF1BP1
Synonyms:	C6orf107; dj349A12.1; ICBP90
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017754
ORF Size:	4320 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221458).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017754.3
RefSeq Size:	9658 bp
RefSeq ORF:	4323 bp
Locus ID:	54887
UniProt ID:	Q6BDS2
Cytogenetics:	6p21.31
MW:	159.3 kDa
Gene Summary:	May act as a negative regulator of cell growth.[UniProtKB/Swiss-Prot Function]



[View online »](#)