

Product datasheet for **RC221469**

Myelin Basic Protein (MBP) (NM_001025101) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Myelin Basic Protein (MBP) (NM_001025101) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: Myelin Basic Protein
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC221469 representing NM_001025101
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAAACCACGCAGGCAAACGAGAATTAATGCCGAGAAGGCCAGTACGAATAGTGAACTAACAGAG
GAGAATCTGAAAAAAGAGAAACCTGGGTGAACTTTCACGGACAACCTCAGAGGACAACGAAGTGTTCGG
AGAGGCAGATGCGAACCAAGAACAATGGGACCTCCTCTCAGGACACAGCGGTGACTGACTCCAAGCGCACA
GCGGACCCGAAGAATGCCTGGCAGGATGCCACCCAGCTGACCCAGGGAGCCGCCCACTTGATCCGCC
TCTTTTCCGAGATGCCCGGGGAGGGAGGACAACACCTTCAAAGACAGGCCCTCTGAGTCCGACGAGCT
CCAGACCATCCAAGAAGACAGTGCAGCCACCTCCGAGAGCCTGGATGTGATGGCGTCACAGAAGAGACCC
TCCCAGAGGCACGGATCCAAGTACCTGGCCACAGCAAGTACCATGGACCATGCCAGGCATGGCTTCTCC
CAAGGCACAGAGACACGGGCATCCTTGACTCCATCGGGCGCTTCTTTGGCGGTGACAGGGGTGCGCCCAA
GCGGGGCTCTGGCAAGGACTCACACCACCCGGCAAGAAGTGTCACTACGGCTCCCTGCCCCAGAAGTCA
CACGGCCGACCCAAGATGAAAACCCCGTAGTCCACTTCTTCAAGAACATTGTGACGCCTCGCACACCAC
CCCCGTCGCAGGAAAGGGGAGAGGACTGTCCCTGAGCAGATTTAGCTGGGGGGCCGAAGCCAGAGACC
AGGATTTGGCTACGGAGGCAGAGCGTCCGACTATAAATCGGCTCACAAGGGATTCAAGGGAGTCGATGCC
CAGGGCAGCTTTCAAAATTTTAAGCTGGGAGGAAGAGATAGTCGCTCTGGATCACCCATGGCTAGAC
GC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC221469 representing NM_001025101
Red=Cloning site Green=Tags(s)

MGNHAGKRELNAEKASTNSETNRGESEKKRNLGELSRRTTSEDNEVFGEADANQNGTSSQDTAVTDSKRT
 ADPKNAWQDAHPADPGSRPHLIRLF SRDAPGREDNTFKDRPSEDELQTIQEDSAATSESLDVMASQKRP
 SQRHGSKYLATASTMDHARHGFLPRHRDTGILDSIGRFGGDRGAPKRGSGKDSHHPARTAHYGLPQKS
 HGRTQDENPVVHFFKNIIVTPRTPPPSQGKGRGLSLSRFSWGAEGQRPFGYGGGRASDYKSAHKGFKGVDA
 QGTL SKIFKLGGDRSRSGSPMARR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8050_h02.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001025101


ORF Size: 912 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001025101.2
RefSeq Size:	2794 bp
RefSeq ORF:	915 bp
Locus ID:	4155
UniProt ID:	P02686
Cytogenetics:	18q23
MW:	32.9 kDa
Gene Summary:	<p>The protein encoded by the classic MBP gene is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the nervous system. However, MBP-related transcripts are also present in the bone marrow and the immune system. These mRNAs arise from the long MBP gene (otherwise called "Golli-MBP") that contains 3 additional exons located upstream of the classic MBP exons. Alternative splicing from the Golli and the MBP transcription start sites gives rise to 2 sets of MBP-related transcripts and gene products. The Golli mRNAs contain 3 exons unique to Golli-MBP, spliced in-frame to 1 or more MBP exons. They encode hybrid proteins that have N-terminal Golli aa sequence linked to MBP aa sequence. The second family of transcripts contain only MBP exons and produce the well characterized myelin basic proteins. This complex gene structure is conserved among species suggesting that the MBP transcription unit is an integral part of the Golli transcription unit and that this arrangement is important for the function and/or regulation of these genes. [provided by RefSeq, Jul 2008]</p>


Product images:



Circular map for RC221469



Western blot validation of overexpression lysate (Cat# [LY422588]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC221469 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified MBP protein (Cat# [TP321469]). The protein was produced from HEK293T cells transfected with MBP cDNA clone (Cat# RC221469) using MegaTran 2.0 (Cat# [TT210002]).