

Product datasheet for **RC221513**

Urotensin II (UTS2) (NM_021995) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Urotensin II (UTS2) (NM_021995) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Urotensin II
Synonyms:	PRO1068; U-II; UCN2; UII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC221513 representing NM_021995 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAACCAACGTATTTTCATCTTATGCTCTGCGTCACTTCTGCTCGGACTCATAAATCCACGTCTCTTT
GCTTTGGCCACTTCAACTCATATCCAAGCCTTCCTTTAATTCATGATTTATTGCTGGAAATATCCTTTCA
ACTCTCAGCACCTCATGAAGACGCGCTTAACCCGGAGGAGCTAGAAAGAGCTTCCTTCTACAGATA
CTGCCAGAGATGCTGGGTGCAGAAAGAGGGATATTCTCAGGAAAGCAGACTCAAGTACCAACATTTTAA
ACCCAAGAGGAAATTTGAGAAAGTTTCAGGATTTCTCTGGACAAGATCCTAACATTTTACTGAGTCATCT
TTTGGCCAGAATCTGAAACCATACAAGAAACGTGAGACTCCTGATTGCTTCTGGAAATACTGTGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC221513 representing NM_021995 Red=Cloning site Green=Tags(s)
-------------------	--

METNVFHLMLCVTSARTHKSTSLCFGHFNSYPSLPLIHDLLLEISFQLSAPHEDARLTPEELERASLLQI
LPEMLGAERGDILRKADSSSTNIFNPRGNLRKFQDFSGQDPNILLSHLLARIWKPKYKRETPDCFWKYCV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:	https://cdn.origene.com/chromatograms/mk8042_g05.zip
Restriction Sites:	Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_021995

ORF Size: 417 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_021995.2
RefSeq Size:	652 bp
RefSeq ORF:	420 bp
Locus ID:	10911
UniProt ID:	O95399
Cytogenetics:	1p36.23
Protein Families:	Secreted Protein
MW:	14 kDa
Gene Summary:	This gene encodes a mature peptide that is an active cyclic heptapeptide absolutely conserved from lamprey to human. The active peptide acts as a vasoconstrictor and is expressed only in brain tissue. Despite the gene family name similarity, this gene is not homologous to urocortin, a member of the sauvagine/corticotropin-releasing factor/urotensin I family. Most of the proprotein is cleaved to make the mature peptide. Transcript variants encoding different preproprotein isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC221513