

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for RC221522L1V

MRP3 (ABCC3) (NM_003786) Human Tagged ORF Clone Lentiviral Particle

Product data:

| Product Type: | Lentiviral Particles |
|------------------------------|---|
| Product Name: | MRP3 (ABCC3) (NM_003786) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MRP3 |
| Synonyms: | ABC31; cMOAT2; EST90757; MLP2; MOAT-D; MRP3 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_003786 |
| ORF Size: | 4581 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC221522). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u> |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | <u>NM 003786.2</u> |
| RefSeq Size: | 5176 bp |
| RefSeq ORF: | 4584 bp |
| Locus ID: | 8714 |
| UniProt ID: | <u>015438</u> |
| Cytogenetics: | 17q21.33 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | ABC transporters |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

| | MRP3 (ABCC3) (NM_003786) Human Tagged ORF Clone Lentiviral Particle – RC221522L1V |
|---------------|--|
| MW: | 169.2 kDa |
| Gene Summary: | The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MRP subfamily which is involved in multi-drug resistance. The specific function of this protein has not yet been determined; however, this protein may play a role in the transport of biliary and intestinal excretion of organic anions. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008] |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US