

Product datasheet for RC221533L4

OTOL1 (NM_001080440) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OTOL1 (NM_001080440) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	OTOL1
Synonyms:	C1QTNF15; C1QTNF16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221533).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001080440
ORF Size:	1431 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001080440.1 , NP_001073909.1
RefSeq Size:	1434 bp
RefSeq ORF:	1434 bp
Locus ID:	131149
UniProt ID:	A6NHN0
Cytogenetics:	3q26.1
Protein Families:	Transmembrane
MW:	49.2 kDa
Gene Summary:	This gene encodes a secreted glycoprotein with a C-terminal complement Cq1-like globular domain that belongs to the C1q/tumor necrosis factor-related protein (CTRP) family. The encoded protein is expressed in the inner ear and forms a multimeric complex called the otoconia, together with cerebellin-1 and otoconin-90, as part of the otoconial membrane. It contains extensive posttranslational modifications including hydroxylated prolines and glycosylated lysines. Naturally occurring mutations in this gene are associated with abnormal otoconia formation and balance deficits resulting from vestibular dysfunction. [provided by RefSeq, Jul 2017]

Product images:



Circular map for RC221533L4