

Product datasheet for **RC221553L4V**

AGRP (NM_007316) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AGRP (NM_007316) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AGRP
Synonyms:	agouti related protein; AGRT; Agrt; ART; Art; ASIP2; MGC118963; orexigenic neuropeptide
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007316
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221553).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007316.1 , NP_015531.1
RefSeq Size:	783 bp
RefSeq ORF:	398 bp
Locus ID:	181
Cytogenetics:	16q22.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Adipocytokine signaling pathway
MW:	14.4 kDa



[View online »](#)

Gene Summary:

This gene encodes an antagonist of the melanocortin-3 and melanocortin-4 receptor. It appears to regulate hypothalamic control of feeding behavior via melanocortin receptor and/or intracellular calcium regulation, and thus plays a role in weight homeostasis. Mutations in this gene have been associated with late on-set obesity. [provided by RefSeq, Dec 2009]