

## Product datasheet for **RC221571L3V**

### Ferredoxin Reductase (FDXR) (NM\_004110) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ferredoxin Reductase (FDXR) (NM_004110) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Ferredoxin Reductase
Synonyms:	ADR; ADXR; ANOA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004110
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221571).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004110.2</a>
RefSeq Size:	1838 bp
RefSeq ORF:	1494 bp
Locus ID:	2232
UniProt ID:	<a href="#">P22570</a>
Cytogenetics:	17q25.1
Protein Families:	Druggable Genome
MW:	54.48 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a mitochondrial flavoprotein that initiates electron transport for cytochromes P450 receiving electrons from NADPH. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Apr 2012]