

Product datasheet for RC221573L3

SRP9 (NM_003133) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: SRP9 (NM_003133) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: SRP9

Synonyms: ALURBP

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC221573).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_003133

ORF Size: 258 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



SRP9 (NM_003133) Human Tagged Lenti ORF Clone - RC221573L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 003133.1</u>

RefSeq Size: 1524 bp
RefSeq ORF: 261 bp
Locus ID: 6726
UniProt ID: P49458

Cytogenetics: 1q42.12

Protein Families: Druggable Genome

Protein Pathways: Protein export

MW: 10.1 kDa


Gene Summary: Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the

rough endoplasmic reticulum membrane. SRP9 together with SRP14 and the Alu portion of the SRP RNA, constitutes the elongation arrest domain of SRP. The complex of SRP9 and

SRP14 is required for SRP RNA binding.[UniProtKB/Swiss-Prot Function]



Product images:



Circular map for RC221573L3