

## Product datasheet for **RC221573L4V**

### SRP9 (NM\_003133) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SRP9 (NM_003133) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SRP9
Synonyms:	ALURBP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003133
ORF Size:	258 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221573).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003133.1</a>
RefSeq Size:	1524 bp
RefSeq ORF:	261 bp
Locus ID:	6726
UniProt ID:	<a href="#">P49458</a>
Cytogenetics:	1q42.12
Protein Families:	Druggable Genome
Protein Pathways:	Protein export



[View online »](#)

**MW:** 10.1 kDa

**Gene Summary:** Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the rough endoplasmic reticulum membrane. SRP9 together with SRP14 and the Alu portion of the SRP RNA, constitutes the elongation arrest domain of SRP. The complex of SRP9 and SRP14 is required for SRP RNA binding.[UniProtKB/Swiss-Prot Function]