

Product datasheet for **RC221580L1V**

CED6 (GULP1) (NM_016315) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CED6 (GULP1) (NM_016315) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CED6
Synonyms:	CED-6; CED6; GULP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016315
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221580).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016315.2
RefSeq Size:	3464 bp
RefSeq ORF:	915 bp
Locus ID:	51454
UniProt ID:	Q9UBP9
Cytogenetics:	2q32.1-q32.2
Domains:	PID
Protein Families:	Druggable Genome



[View online »](#)

MW: 34.3 kDa

Gene Summary: The protein encoded by this gene is an adapter protein necessary for the engulfment of apoptotic cells by phagocytes. Several transcript variants, some protein coding and some thought not to be protein coding, have been found for this gene. [provided by RefSeq, Nov 2011]