

Product datasheet for **RC221641L4V**

HDLBP (NM_005336) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HDLBP (NM_005336) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HDLBP |
| Synonyms: | HBP; PRO2900; VGL |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_005336 |
| ORF Size: | 3804 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC221641). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_005336.2 |
| RefSeq Size: | 5032 bp |
| RefSeq ORF: | 3807 bp |
| Locus ID: | 3069 |
| UniProt ID: | Q00341 |
| Cytogenetics: | 2q37.3 |
| Domains: | KH |
| MW: | 141.3 kDa |


[View online »](#)

Gene Summary:

The protein encoded by this gene binds high density lipoprotein (HDL) and may function to regulate excess cholesterol levels in cells. The encoded protein also binds RNA and can induce heterochromatin formation. [provided by RefSeq, Mar 2016]