

Product datasheet for **RC221662**

RAI16 (FAM160B2) (NM_022749) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAI16 (FAM160B2) (NM_022749) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RAI16
Synonyms:	RAI16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC221662 representing NM_022749
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGCTGAGCCGGCTCGGGCGCTGCTGCAGGAAGCCGTGGGGCGCGGAGCCCAGCATTGACCTGCTGC
 AGGCCITTCGTGGAGCACTGGAAGGGCATCACGCACTACTACATCGAGAGCACAGATGAAAGCACCCCGC
 CAAGAAGACAGACATTCCCTGGCGGCTGAAGCAGATGCTGGATATCCTGGTGTATGAAGAGCAGCAGCAG
 GCGGCCGCGGGTGAAGCAGGGCCCTGCCTGGAGTACCTGCTGCAGCACAAGATCCTGGAGACTCTCTGCA
 CGCTGGGCAAGGCCGAGTACCCCCAGGCATGCGGCAGCAGGTGTTCCAGTTCTTCAGCAAGGTTCTGGC
 GCAGGTGCAGCACCCCTGCTGCATTACCTCAGCGTCCACAGGCCTGTGCAGAACTCCTCCGACTTGGT
 GGGACTGCTTCCGGATCCGTTACAGAAAAGGAGGAGGTGCAGTTCACCACCGTCTCTGCTCCAAGATCC
 AGCAGGACCCAGAGCTGCTCGCTACATCCTGGAAGGTAAGAAAGATTGTAGGTAGGAAGAAAGCATGCGG
 AGAACCCACTGCCCTGCCTAAGGACACAACCAGCCACGGGGACAAGGACTGCTCCACGATGGTGTCTCT
 GCCAGGCCCCAGCTGGACGGGGAGTCTGTGGGGCCAGGCCCTGAACAGCCACATGCCTGCTGAGACCG
 AGGAGCTGGACGGTGGGACCACAGAGAGCAACCTGATTACCTCCCTGCTTGGGCTGTGCCAGAGCAAGAA
 GAGTCGGGTGGCCTTGAAGGCCAGGAGAACCTGCTGCTCCTGGTGAAGATGGCCTCCCCAGCAGCTGCC
 ACCTACCTGGTACAGAGCAGCGCTGCTGCCCTGCGATCGTCCGGCACCTTTGCCAGTTGTACCGGTCCA
 TGCCTGTCTTCTGGACCCCGCAGACATTGCCACCTTAGAGGGCATCAGCTGGAGGTTACCCAGTGCCCC
 GTCTGATGAGGCTTCTTCCCTGGCAAGGAGGCCTTGGCTGCCTTCTTGGGCTGGTTTGATTACTGCGAC
 CACCTCATCACAGAGGCACACAGGTGGTTGCGGACGCCTTGGCAAGGCTGTGGCTGAGAACTTCTTCG
 TGGAGACCCTGCAGCCCCAGCTCCTGCACGTGTCCGAGCAGAGCATCTTGACCTCCACCGCCCTCCTCAC
 AGCCATGCTGCGCCAGCTTCGCTCCCTGCGCTGCTGCGGGAGGCGGTGGCTTTCTCTCTGGGCACAGAC
 CGGCAGCCTGAAGCCCCGGGGACAACCCACACCCCTGTATGCTCATCTCATCGGGCATTGTGACCACC
 TCTCTGATGAGATCAGCATCACCACTCCGGCTGTTTGAGGAGCTGCTGCAGAAGCCCCACGAGGGGAT
 CATCCACAGCCTGGTCTGCGCAACCTTGAGGGCCGCCCTTACGTGGCCTGGGGCTACCCAGAGCCTGAG
 AGCTATGAGGACACCCTAGACCTGGAGGAAGACCCCTACTTCACCGACAGCTTCTGGATTCCGGCTTTC
 AAACCCCGCAAAGCCTCGCTAGCTCCTGCTACAGTTACGATGGCAAAACAGCAGTGACCGAGATCGT
 CAACAGTTTCTGTGCCTGGTCCCCGAGGAAGCCAAGACCTCTGCCTTCTGGAGGAGACAGGCTATGAC
 ACATACGTCCACGATGCTTATGGCTGTTCAGGAGTGCAGCTCCCGCTGCGCTCCTGGGGCTGGCCTC
 TGACCCCCACACCTTTGGACCCCATGAGCCCGAGCGACCTTTCTTCGAGGGCCACTTCTCCGAGTGCT
 GTTTGACCGCATGTCCCGGATTCTGGATCAGCCATACAGCCTGAACCTGCAGGTGACCTCGGTCTGTCC
 CGGCTTGGCCTCTTCCCCACCCCATATTATGAGTACCTGCTGGATCCGTACATCAGCCTGGCCCCCG
 GCTGCAGGAGCCTATTCTCCGTGTTGGTGAAGGTGATCGGGGACTTGATGCAGAGAATCCAGAGGGTACC
 CCAGTTCCCAGGCAAGCTGCTCCTGGTGCAGCAGTTCAGGGGCCAGGCTCCTGGGGAGCAGCTGGAC
 CACCAGACCCTCCTCCAGGGCGTGGTGGTGGTGGAGGAGTTCTGCAAGGAGCTGGCTGCCATTGCCTTCG
 TCAAGTTTCCCCACATGATCTCGCCAGAAGCTCTCCCAGCCCCGAAGGGCAGGT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC221662 representing NM_022749
 Red=Cloning site Green=Tags(s)


MLSRLGALLQEAVGAREPSIDLLQAFVEHWKGITHYYIESTDESTPAKKTIPWRLKQMLDILVYEEQQQ
 AAAGEAGPCLEYLLQHKILETLCTLGKAEYPPGMRQQVFQFFSKVLAQVQHPLLHYLSVHRPVQKLLRLG
 GTASGSVTEKEEVQFTTVLCSKIQQDPELLAYILEGKKIVGRKKACCEPTALPKDTTSHGDKDCSHDGAP
 ARPQLDGESCGAALNSHMPAETEELDGGTTESNLITSLGLCQSKSRVALKAQENLLLLVSMASPA
 TYLVQSSACCPAIVRHLCQLYRSMVPVFLDPADIATLEGISWRLPSAPSDEASFPGKEALAAFLGWFDYCD
 HLITEAHTVVADALAKAVAENFFVETLQPQLLVHVEQSILTSTALLTAMLRQLRSPALLREAVAFLLGTD
 RQPEAPGDNPHLYAHLIGHCDHLSDEISITTLRLFEELLQKPHEGIIHSLVLRNLEGRPYVAWGSPEPE
 SYEDTLDLEEDPYFTDSFLDSGFQTPAKPRLAPATSYDGKTAVTEIVNSFLCLVPEEAKTSAFLEETGYD
 TYVHDAYGLFQECSSRVASWGWPLTPTLPDHEPERPFFEGHFLRVLFDRMSRILDQPYSLNLQVTSVLS
 RLALFPHPIHEYLDPYISLAPGCRSLFSVLVRVIGDLMQRIQRPQFPKLLLVRKQLTGQAPGEQLD
 HQTLLQGVVVLEEFCKELAAIAFVKFPPHDPQNVSPAPEGQV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8040_a01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022749


ORF Size: 2229 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_022749.7</u>
RefSeq Size:	3759 bp
RefSeq ORF:	2232 bp
Locus ID:	64760
UniProt ID:	<u>Q86V87</u>
Cytogenetics:	8p21.3
MW:	82.2 kDa
Gene Summary:	Able to activate MAPK/ERK and TGFB signaling pathways (PubMed:22971576). May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation (By similarity). May play a role in cell proliferation (PubMed:22971576). [UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC221662



Western blot validation of overexpression lysate (Cat# [LY411587]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC221662 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified FAM160B2 protein (Cat# [TP321662]). The protein was produced from HEK293T cells transfected with FAM160B2 cDNA clone (Cat# RC221662) using MegaTran 2.0 (Cat# [TT210002]).