


Product datasheet for RC221663L4

EPB41L4A (NM_022140) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EPB41L4A (NM_022140) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	EPB41L4A
Synonyms:	EPB41L4; NBL4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221663).
Restriction Sites:	SgfI-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_022140
ORF Size:	2058 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022140.3 , NP_071423.3
RefSeq Size:	4707 bp
RefSeq ORF:	2061 bp
Locus ID:	64097
UniProt ID:	Q9HCS5
Cytogenetics:	5q22.1-q22.2
Domains:	B41
MW:	78.9 kDa
Gene Summary:	The protein encoded by this gene is a member of the band 4.1 protein superfamily. Members of this superfamily are thought to play an important role in regulating interactions between the cytoskeleton and plasma membrane, and contain an amino terminal conserved domain that binds glycophorin C. This gene product is thought to be involved in the beta-catenin signaling pathway. [provided by RefSeq, Dec 2016]

Product images:



Circular map for RC221663L4