

Product datasheet for **RC221679A1V**

Human CGB2 (NM_033378) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CGB2 (NM_033378) AAV Particle
Tag: Myc-DDK
Symbol: CGB2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC221679 representing NM_033378
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCAAAGGGGCTGCTGCTGTTGCTGCTGCTGAGCATGGGCGGGACATGGGCATCCAAGGAGCCGCTTC
GGCCACGGTGCCGCCCATCAATGCCACCCTGGCTGTGGAGAAGGAGGGCTGCCCGTGTGCATCACCGT
CAACACCACCATCTGTGCCGGCTACTGCCCCACCATGACCCGCGTGTGCAGGGGGTCTGCCGGCCCTG
CCTCAGGTGGTGTGCAACTACCGCGATGTGCGCTTCGAGTCCATCCGGTCCCTGGCTGCCCGCGCGGCG
TGAACCCCGTGGTCTCCTACGCCGTGGCTCTCAGCTGTCAATGTGCACTCTGCCCGCGCAGCACCCTGA
CTGCGGGGGTCCAAGGACCACCCCTTGACCTGTGATGACCCCGCTTCCAGGCCCTCTTCTCTCAAAG
GCCCTCCCCCAGCCTTCCAAGCCATCCCGACTCCCGGGGCCCTCAGACACCCGATCTCCACAA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC221679 representing NM_033378
Red=Cloning site Green=Tags(s)

MSKGLLLLLLLSMGGTWASKEPLRPRCPINATLAVEKEGCPVCITVNTTICAGYCPTMTRVLQGVLPAL
PQVVCNYRDVRFESIRLPGCPRGVNPVSYAVALSCQCALCRRSTTDCGGPKDHPLTCDPRFQASSSSK
APPPSLPSPRLPGPSDTPILPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_033378



[View online »](#)

ORF Size:	489 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033378.1</u>
RefSeq Size:	732 bp
RefSeq ORF:	492 bp
Locus ID:	114336
UniProt ID:	<u>Q6NT52</u>
Cytogenetics:	19q13.33
MW:	17.2 kDa