

Product datasheet for **RC221687**

BCAS4 (NM_001010974) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: BCAS4 (NM_001010974) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: BCAS4
Synonyms: CNOL
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC221687 representing NM_001010974
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGCGGACCGGGGGCGGGCTCCGAGGCCCGGGCGCAACCACGGGCTCCCAGGCAGCCTCCGCCAGC
CGGACCCCGTCGCCCTCCTGATGCTGCTCGTGGACGCTGATCAGCCGGAGCCCATGCGCAGCGGGGCGCG
CGAGCTCGCGCTTCTGACCCCGAGCCTGGGGCCGAGGCGAAGGAGGTGGAGGAGACCATCGAGGGC
ATGCTCCTCAGGCTGGAAGAGTTTTGCAGCCTGGCTGACCTGATCAGGAGTGATACTTCACAGATCTGG
AGGAAAACATCCCAGTCCTTAAGGCCAACTGACAGAAATGCGTGGCATCTATGCCAAAGTGGACCGGCT
AGAGAAGTCACCTGCACCGGTGCCCGTACGTACGAGCTGCCACACTGTATAGGACGGAGGACTATTTT
CCTGTGGACGCCGGGAAGCACAGCACACCACCCCGCACCTGCCCTCGGCCTTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC221687 representing NM_001010974
Red=Cloning site **Green**=Tags(s)

MQRTGGGAPRPNHGLPGSLRQDPVALLMLLVDADQPEPMRSGARELALFLTPEPGAEAKEVEETIEG
MLLRLEEFCSLADLIRSSTSQILEENIPVLKAKLTEMRGIYAKVDRLEKSPAPVPVTYELPTLYRTEDYF
PVDAGEAQHPRTCPRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1461_a08.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001010974

ORF Size: 474 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001010974.2](#), [NP_001010974.1](#)


RefSeq Size: 1136 bp

RefSeq ORF: 477 bp


Locus ID: 55653

UniProt ID: Q8TDM0
 Cytogenetics: 20q13.13
 MW: 17.3 kDa

Product images:



Circular map for RC221687



Western blot validation of overexpression lysate (Cat# [LY423248]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC221687 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).