

Product datasheet for **RC221699L1V**

CD105 (ENG) (NM_000118) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	CD105 (ENG) (NM_000118) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD105
Synonyms:	END; HHT1; ORW1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000118
ORF Size:	1884 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000118.1
RefSeq Size:	3142 bp
RefSeq ORF:	1878 bp
Locus ID:	2022
UniProt ID:	P17813
Cytogenetics:	9q34.11
Domains:	zona_pellucida
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane


[View online »](#)

MW: 67.78 kDa

Gene Summary: This gene encodes a homodimeric transmembrane protein which is a major glycoprotein of the vascular endothelium. This protein is a component of the transforming growth factor beta receptor complex and it binds to the beta1 and beta3 peptides with high affinity. Mutations in this gene cause hereditary hemorrhagic telangiectasia, also known as Osler-Rendu-Weber syndrome 1, an autosomal dominant multisystemic vascular dysplasia. This gene may also be involved in preeclampsia and several types of cancer. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2013]

Product images:



[RC221699L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC221699L1V particle to overexpress human ENG-Myc-DDK fusion protein.