

## Product datasheet for **RC221739L3V**

### TPD52 (NM\_001025252) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TPD52 (NM_001025252) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TPD52
Synonyms:	D52; hD52; N8L; PC-1; PrLZ
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025252
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025252.1</a>
RefSeq Size:	4250 bp
RefSeq ORF:	675 bp
Locus ID:	7163
UniProt ID:	<a href="#">P55327</a>
Cytogenetics:	8q21.13
MW:	24.1 kDa



[View online »](#)