

Product datasheet for **RC221750A1V**

Human OR6Q1 (NM_001005186) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR6Q1 (NM_001005186) AAV Particle
Tag:	Myc-DDK
Symbol:	OR6Q1
Synonyms:	OR11-226
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221750 representing NM_001005186 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAACCATATACCAAAACTGGACCCAGGTAAGTGAATTTGTCATGATGGGCTTTGCTGGCATCCATG
AAGCACACCTCCTCTTCTTCACTCTTCCTCACCATGTACCTGTTACCTTGGTGGAGAATTTGGCCAT
CATTTTAGTGGTGGGTTTGGACCACCGACTACGGAGACCCATGTATTTCTTCTGACACACTTGTCTGC
CTTGAAATCTGGTACACTTCTGTTACAGTGCCCAAGATGCTGGCTGGTTTATTGGGGTGGATGGTGGCA
AGAATATCTCTTATGCTGGTTGCCTATCCAGCTCTTCATCTTCACCTTTCTTGGGGCACTGAGTGT
CCTACTGGCTGCCATGGCCTATGATCGTTATGTGGCCATTTGTATGCCTCTCCACTATGGGGCTTTTGTG
TCCTGGGGCACCTGCATCCGTCTGGCAGCTGCCTGTTGGCTGGTAGGTTTCCTCACACCCATCTTGCCAA
TCTACCTCTTGTCTCAGCTAACATTTTGTGGCCAAATGTCATTGACCATTTCTCCTGTGATGCCTCACC
CTTGCTAGCCTTGTCTGCTCAGATGTCACTTGAAGGAGACTGTGGATTTCTGGTGTCTCTGGCTGTG
CTACTGGCCTCCTCTATGGTCATTGCTGTGCTATGGCAACATCGTCTGGACACTGCTGCACATCCGCT
CAGCTGCTGAGCGCTGGAAGGCCTTCTACCTGTGCAGCTCACCTGACTGTGGTGAAGCCTCTTCTATGG
CACTCTTTTCTTATGTATGTCCAGACCAAGGTGACCTCCTCCATCAACTCAACAAGGTGGTATCTGTC
TTCTACTCTGTTGTACGCCCATGCTCAATCCTCTCATCTACAGTCTTAGGAACAAGGAAGTGAAGGGAG
CTCTGGGTCGAGTCTTTTCTCTCAACTTTTGAAGGGACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC221750 representing NM_001005186 Red=Cloning site Green=Tags(s) MQPYTKNWTQVTEFVMMGFAGIHEAHLFFILFLTMYLFTLVENLAIILVVGLDHRLRRPMYFFLTHLSC LEIWYTSVTVPKMLAGFIGVDGGKNISYAGCLSQLFIFTFLGATECFLLAAMAYDRYVAICMPLHYGAFV SWGTCIRLAAACWLVGFLTPILPIYLLSQLTFCGPNVIDHFSCDASPLLALSCSDVTWKETVDFLVSLAV LLASSMVIASVYGNIVWTLHLIRSAERWKAFSTCAHLTVVSLFYGTLLFFMYVQTKVTSSINFNKVSV FYSVVTPLNPLIYSLRNKEVGALGRVFSLNFWKGQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001005186
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005186.1, NP_001005186.1</u>
RefSeq Size:	954 bp
RefSeq ORF:	954 bp
Locus ID:	219952
UniProt ID:	<u>Q8NGQ2</u>
Cytogenetics:	11q12.1
MW:	35.6 kDa