

Product datasheet for **RC221774L4V**

ART5 (NM_001079536) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ART5 (NM_001079536) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ART5
Synonyms:	ARTC5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001079536
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221774).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001079536.1 , NP_001073004.1
RefSeq Size:	1249 bp
RefSeq ORF:	876 bp
Locus ID:	116969
UniProt ID:	Q96L15
Cytogenetics:	11p15.4
Protein Families:	Secreted Protein
MW:	32.2 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the ARG-specific ADP-ribosyltransferase family. Proteins in this family regulate the function of target proteins by attaching ADP-ribose to specific amino acid residues in their target proteins. The mouse homolog lacks a glycosylphosphatidylinositol-anchor signal sequence and is predicted to be a secretory enzyme. Several transcripts encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]